

Product Description

Name of the Device (s): AP | Instrument Protector

Product (group) code (s): BM3300258*, BM3300259*

Description:

The Instrument Protector is intended to be used as a protecting device to prevent sharp medical instruments from damaging the sterile barrier or hurting the end-user. The Instrument Protector is only intended to be used on small, critical surgical instruments.

The sterilization method which suits this product is STEAM. The Instrument Protector is for single use. The Instrument Protector is conform the standard ISO 11607-1:2020.



Place, date (dd.mm.yyyy) of issue and version:

Zaandam, 24.11.2022, v5