

Technical Data Sheet

Manufacturers Name:	Blue Medical B.V.
Manufacturers Address:	Sluispolderweg 79, 1505HJ, Zaandam, The Netherlands
SRN (Single Registration Number):	NL-MF-000010818
Name of the Device (s):	AP Instrument Protector
Product group code (s):	BM3300258*, BM3300259*
Classification:	CLASS I Medical Device
Manufacturer product category:	AP
Manufacturer product group:	AP Instrument Protector
Market product category:	Ancillary products
Product material:	Medical grade paper + PET foil
Product color:	White
Product size:	BM33002580: 51mm x 127mm BM33002585: 89mm x 168mm BM33002590: 51mm x 51mm BM33002595: 56mm x 86mm
Applicable standards:	ISO 11607-1:2020
Suitable processes:	STEAM
Standard process values:	134°C – 3,5 min., 134°C – 18 min.
Instructions for use:	YES, eIFU www.blue-medical.nl/eifu
HS / Customs code:	3822.9000
Shelf life:	5 years
Lead free:	YES
Latex free:	YES

DEHP / CMR free:

YES

Storage conditions:

Humidity limitation 30% - 70% / Temperature limit +10°C - +30°C / Keep dry / Keep away from sunlight

Product symbol markings:

Do not re-use / Do not use if package is damaged / CE / Catalogue number (REF) / Batch code (LOT) / Consult instructions for use / Date of manufacture / Humidity limitation / Keep away from sunlight / Manufacturer / Medical Device / STEAM / Temperature limit / Use-by date / Keep dry

Box content:

100 pieces

Additional note:

-

Articles valid under this TDS:

Article number	UDI-DI (GTIN)	Description	Box content
BM33002580	8720174715631	AP Instrument Protector small – 100 pcs	100 pcs.
BM33002585	8720174715648	AP Instrument Protector medium – 100 pcs	100 pcs.
BM33002590	8720174715655	AP Instrument Protector mini – 100 pcs	100 pcs.
BM33002595	8720174715662	AP Instrument Protector 56mm x 86mm – 100 pcs	100 pcs.

This Technical Data Sheet is issued under the sole responsibility of Blue Medical B.V.. We hereby declare that the product(s) specified above meet the indicated standards. All supporting documentation is retained at the premises of the manufacturer.

Signature:



Place, date (dd.mm.yyyy) of issue and version:

Mr. Alex Welleweerd
Sales Director

Zaandam, 24.02.2023, v8