

Product Description

Name of the Device (s): AP | Seal Check

Product (group) code (s): BM3300256*

Description:

The Seal Check is intended to be used as verification tool to verify the sealing quality performed by a sealing machine. By inserting the Seal Check into the pouch/reel and by performing a seal, a clear indication of the seal quality will appear as result of the contrast pattern of the Seal Check. The width of the Seal Check is 175mm which ensures a complete coverage of the sealing wheel. All relevant information such as sealing machine number, location, date and results can be noted on the Seal Check.

The Seal Check with reference BM33002560 can be used on film/paper seal pouches and reels. The Seal Check NW Polyolefin fabric with reference BM33002561 can be used on film/polyolefin fabric (Tyvek®) pouches and reels. The Seal Check is not to be sterilized and is for single use. The Seal Check supports the standard ISO 11607-2:2020.



Place, date (dd.mm.yyyy) of issue and version:

Zaandam, 19.04.2022, v3