

Technical Data Sheet

Manufacturers Name:	Blue Medical B.V.
Manufacturers Address:	Sluispolderweg 79, 1505HJ, Zaandam, The Netherlands
SRN (Single Registration Number):	NL-MF-000010818
Name of the Device (s):	AP Seal Check
Product group code (s):	BM3300256*
Classification:	Non-Medical Device
Manufacturer product category:	AP
Manufacturer product group:	AP Seal Check
Market product category:	Ancillary products
Product material:	BM33002560: Medical grade paper (70g/m ²) BM33002561: Tyvek® (75g/m ²)
Product color:	N/A
Product size:	175mm x 75mm
Applicable standards:	N/A
Instructions for use:	YES, eIFU www.blue-medical.nl/eifu
HS / Customs code:	3822.9000
Shelf life:	5 years
Lead free:	YES
Latex free:	YES
DEHP / CMR free:	YES
Storage conditions:	Humidity limitation 30% - 70% / Temperature limit +10°C - +30°C / Keep dry / Keep away from sunlight

Product symbol markings:

Do not re-use / Do not use if package is damaged / Catalogue number (REF) / Batch code (LOT) / Consult instructions for use / Date of manufacture / Humidity limitation / Keep away from sunlight / Manufacturer / Temperature limit / Use-by date / Keep dry

Box content:

See article summary below.

Additional note:

-

Articles valid under this TDS:

Article number	UDI-DI (GTIN)	Description	Box content
BM33002560	8720174715563	AP Seal Check 175mm – 250 pcs	250 pcs
BM33002561	8720174715570	AP Seal Check NW Polyolefin fabric 175mm - 250 pcs	250 pcs

This Technical Data Sheet is issued under the sole responsibility of Blue Medical B.V.. We hereby declare that the product(s) specified above meet the indicated standards. All supporting documentation is retained at the premises of the manufacturer.

Signature:



Place, date (dd.mm.yyyy) of issue and version:

Mr. Alex Welleweerd
Sales Director

Zaandam, 17.02.2023, v8