Steri24

Best Prices. Top Service.



Autoclave from
€ 1,399

Sealer from
€ 475

WD from
€ 7,887

ON-SITE SERVICE

3 YEAR WARRANTY

AI-generated image

valid until 31.12.2027

Table of Contents

Preface 3
Instrument Reprocessing 4

Autoclaves

Class B Autoclave Pro 3 L 7
Class B Autoclave Pro 12 L 8
Class B Autoclave Pro 18 L 9
Class B Autoclave Pro 23 L 10
Class B Autoclave Premium 12 L 11
Class B Autoclave Premium 18 L 12
Class B Autoclave Premium 23 L 13
Autoclave Accessories 14

Ultrasonic Cleaners

Steri24 Ultrasonic Cleaners 17

Thermal Disinfector

Miele Thermal Disinfector PWD 8534 19
Miele Thermal Disinfector PWD 8682 20
Miele Accessories 21

Packaging & Sealing Devices

CertoSeal Pro Touch 23
CertoSeal 200 24
CertoSeal 300 25
Sealing Accessories 26
Sterilization Containers 27

Services 28
Contact 29

Steri24 – a brand of CertoClav: A guarantee of quality since 1957

Dear Customers,

a warm welcome to Steri24, your leading provider of medical Class B autoclaves per EN13060! As Managing Director of this long-standing family business in Linz, I personally stand for the highest quality and excellent service.
For over 65 years, CertoClav has stood for reliability, innovation and commitment in medical sterilization technology. Our success is built on close collaboration with you, our valued customers. Our main concern is to offer you the best quality available on the market.

Our Class B autoclaves meet the strict requirements of EN13060 and even exceed them. This allows us to deliver top products that guarantee the highest safety standards.
Our autoclaves are made from the finest materials and undergo thorough quality checks to ensure flawless functionality. You can also rely on our expert team, who are extensively trained in validation and maintenance. Should any questions or issues arise, our excellent on-site service in Germany, Austria and Switzerland is there for you.

Another reason to choose Steri24 is our unbeatable prices. We believe high-quality medical equipment should be affordable. That is why we offer not only top quality, but also fair terms.
As a mark of our confidence in our products, all our autoclaves come with a 2-year warranty. Your satisfaction drives us, and we want you to be completely happy with your purchase.
And that is not all! We also offer free delivery to Germany, Austria and Switzerland. Wherever you are, we get your order to you safely and quickly.

Our commitment to your safety and satisfaction does not end with the purchase. Our 24-hour emergency hotline is always available should you need support. Because at Steri24, you are not just a customer – you are part of our family.
See for yourself the quality, unbeatable prices and outstanding service that Steri24 offers. Let us set the highest standards in medical sterilization together and build a healthier future.

Thank you for your trust.

Yours, Michael Simon, M.Sc.
Managing Director, Steri24



The benefits of Steri24:

Free delivery in DE, AT and CH
for orders of € 1,200 net or more


2–3 year warranty
2 years' warranty as standard; 3 years with a maintenance contract and maintenance carried out per the contract

24-hour emergency hotline
+49 545 263 499 31 - Steri24 service location near you

On-site service
Commissioning, maintenance and validation of your equipment on-site by our service partners

CertoSeal 300

validatable sealer with autocutter



3

Packaging



Miele Thermal Disinfector

Disinfection & drying unit

2

Disinfection & Drying



Cleaning



Steri24 Ultrasonic Cleaners

Available in various sizes

3–22 L

1



AQUAdem 40

Water treatment system



Class B Autoclave Pro or Premium

Available as 3–23 litre units

4 Sterilization



SteriTrace Pro Assistant

Paperless documentation



5 Traceability



Tests & Indicators

For monitoring & control



Label printer

Batch labelling

STERILIZATION
CLASS B AUTOCLAVES



AI-generated image

Class B Autoclave Pro 3 L

€ 1,399

The Space Saver - Lightning fast 11 minutes

- ✓Class B autoclave per EN 13060
- ✓Especially compact for small workspaces
- ✓**Lightning fast**: 7 min for unwrapped instruments, 15 min for wrapped instruments
- ✓Digital SD card documentation & integrated protocol printer
- ✓Color touchscreen for easy operation
- ✓3-litre chamber volume with 2 trays
- ✓No external water connection required



The Steri24 Class B Autoclave Pro Touch 3 L is a compact steam sterilizer for professional instrument reprocessing in the smallest of spaces. It is suitable for wrapped and unwrapped instruments, textiles and hollow items, and meets Class B requirements per EN 13060.

The large color touchscreen makes operation simple and intuitive. In the quick program, unwrapped instruments are ready for use again after just 7 minutes. Integrated SD card and thermal printer documentation ensures complete, compliant logging of sterilization processes.

With a 3-litre chamber volume, 2 trays, an integrated fresh water tank and a compact footprint, the Pro Touch 3 L is ideal for small practices and facilities with limited space. It is especially suited to podiatry and foot care practices, tattoo and piercing studios, beauty and cosmetic studios, and small medical practices.



Scope of delivery:
Steri24 Class B Autoclave Pro 3 L with tray holder
2x quick-connect adapters and 1x hose for manually draining the feed tank and waste water
2x trays, dimensions: 28.5 x 18.5 x 1.5 cm (L x W x H)
Handle for safe tray removal after sterilization
2x spare fuses, 5mm, 10A
Bacteria filter, 1x thermal paper
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Pro 3 L			
Item Number	8500491	**Feed Water Tank Volume**	1 litre (for at least 5 cycles)
Chamber Volume	3 litres	**Sterilization Temperature**	121 °C and 134 °C
Chamber Dimensions (L x W x H)	32 x 19.5 x 7.5 cm	**Outer Dimensions (L x W x H)**	59.5 x 28.5 x 24.5 cm
Documentation	SD card and thermal printer	**Power Connection**	230V AC, 50 Hz
Power	2900 W	**Net Weight**	24 kg

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300557
Thermal paper rolls (10 pcs.)	1300742
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Water distiller	8500482

Class B Autoclave Pro 12 L

€ 1,499

For small & medium-sized practices

- ✓Class B autoclave per EN 13060
- ✓For wrapped and unwrapped instruments
- ✓Stand-alone unit with integrated fresh & waste water tank
- ✓Digital USB documentation & integrated protocol printer
- ✓User-friendly menu navigation
- ✓12-litre chamber volume with 3 trays
- ✓Optional fixed water connection



The Steri24 Class B Autoclave Pro 12 L is a compact steam sterilizer for professional instrument reprocessing in practices and clinics. It is suitable for wrapped and unwrapped instruments, hollow items and textiles, and supports Class B, S and N programs.

Thanks to the integrated fresh and waste water tanks, the unit can be used flexibly as a stand-alone solution. No external water connection is required. Operation is simple and clear, making the autoclave ideal for daily practice use.

With a 12-litre chamber volume, 3 trays, USB documentation, an integrated protocol printer and an optional fixed water connection, the Pro 12 L offers a reliable, compliant solution for small to medium-sized practices. It is especially suited to general medicine, dermatology, ENT, podiatry and foot care practices, and smaller dental and veterinary practices.



Scope of delivery:
Steri24 Class B Autoclave Pro 12 L
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
3x trays, dimensions: 31 x 15 x 1.5 cm (L x W x H)
Tray holder
Handle for safe tray removal after sterilization
Bacteria filter, 1x thermal paper
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Pro 12 L			
Item Number	8500475	**Feed Water Tank Volume**	3 litres (for at least 5 cycles)
Chamber Volume	12 litres	**Waste Water Tank Volume**	3 litres
Chamber Dimensions	Ø 20 cm x 35 cm	**Outer Dimensions (L x W x H)**	64 x 44 x 36 cm
Sterilization Temperature	121°C and 134°C	**Power Connection**	230V AC, 50 Hz
Power	1950 W	**Net Weight**	51 kg
Documentation	USB and thermal printer		

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300557
Thermal paper rolls (10 pcs.)	1300742
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Water distiller	8500482
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562


Class B Autoclave Pro 18 L

€ 1,649

For medium & large practices

- ✓Class B autoclave per EN 13060
- ✓For wrapped and unwrapped instruments
- ✓Suitable for textiles and hollow items
- ✓Stand-alone unit with integrated fresh and waste water tank
- ✓USB documentation and integrated protocol printer
- ✓18-litre chamber volume with 5 stainless steel trays
- ✓Dual drying via vacuum and jacket heating
- ✓Optional fixed water connection





The Steri24 Class B Autoclave Pro 18 L is the perfect steam sterilizer for instrument reprocessing in practices and clinics. The renowned notified body TÜV SÜD carries out the tests and assessments for the conformity assessment and certifies compliance according to uniform assessment standards.

The Steri24 Class B Autoclave Pro stands out above all for its simple operation, durability, safety and digital documentation. Thanks to international production and high volumes, the price is competitive.

The Steri24 Class B Pro is suitable for clinics, dental practices, piercing and tattoo studios, podiatry and foot care practices, cosmetic and eye clinics.

It is a high-quality steam sterilizer that meets strict medical requirements.

Scope of delivery:
Steri24 Class B Autoclave Pro 18 L
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
5x trays, dimensions: 28 x 19 x 1.5 cm (L x W x H)
Tray holder
Handle for safe tray removal after sterilization
Bacteria filter, 1x thermal paper
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Pro 18 L			
Item Number	8500547	**Feed Water Tank Volume**	4 litres (for at least 5 cycles)
Chamber Volume	18 litres	**Waste Water Tank Volume**	4 litres
Chamber Dimensions	Ø 25 cm x 35 cm	**Outer Dimensions (L x W x H)**	64 x 50 x 43 cm
Sterilization Temperature	105 °C to 134 °C	**Power Connection**	230V AC, 50 Hz
Power	1950 W	**Net Weight**	53 kg
Documentation	USB and thermal printer		

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300557
Thermal paper rolls (10 pcs.)	1300742
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Water distiller	8500482
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562

Class B Autoclave Pro 23 L

For large practices & clinics

€ 1,949


- ✓Class B autoclave per EN 13060
- ✓For wrapped and unwrapped instruments
- ✓Suitable for textiles and hollow items
- ✓Stand-alone unit with integrated fresh & waste water tank
- ✓USB documentation and integrated protocol printer
- ✓23-litre chamber volume with 5 trays
- ✓Dual drying via vacuum and jacket heating
- ✓Optional fixed water connection



The Steri24 Class B Autoclave Pro 23 L is a powerful steam sterilizer for professional instrument reprocessing in practices, healthcare centres and clinics. It is suitable for wrapped and unwrapped instruments, textiles and hollow items, and meets Class B requirements per EN 13060.

With a 23-litre chamber volume and 5 stainless steel trays, the autoclave offers extra capacity for facilities with a high volume of instruments. The 7 standard programs, simple operation and dual drying ensure safe, well-dried loads in daily use.

Thanks to the integrated fresh and waste water tanks, the unit can be used flexibly as a stand-alone solution. USB documentation and an integrated protocol printer ensure complete documentation of sterilization processes. The optional fixed water connection provides added flexibility for daily use. Ideal for dental practices, clinics, healthcare centres, surgical practices and veterinary clinics.



Scope of delivery:
Steri24 Class B Autoclave Pro 23 L
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
5x trays, dimensions: 38 x 19 x 2 cm (L x W x H)
Tray holder
Handle for safe tray removal after sterilization
Bacteria filter, 1x thermal paper
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Pro 23 L			
Item Number	8500477	**Feed Water Tank Volume**	4 litres (for at least 5 cycles)
Chamber Volume	23 litres	**Waste Water Tank Volume**	4 litres
Chamber Dimensions	Ø 25 cm x 45 cm	**Outer Dimensions (L x W x H)**	67 x 50 x 44 cm
Sterilization Temperature	105 °C to 134 °C	**Power Connection**	230V AC, 50 Hz
Heating Power	1950 W	**Net Weight**	55 kg
Documentation	USB and thermal printer		

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300557
Thermal paper rolls (10 pcs.)	1300742
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Water distiller	8500482
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562

Class B Autoclave Premium 12 L

€ 2,449

with automatic documentation

- ✓Class B autoclave per EN 13060
- ✓SteriTrace Pro Assistant for digital batch release
- ✓Paperless documentation of sterilization processes
- ✓LAN/WLAN, cloud and practice software connectivity
- ✓Integrated fresh & waste water tank
- ✓Label printer included
- ✓12-litre chamber volume with 3 trays
- ✓Optional fixed water connection



The Steri24 Class B Autoclave Premium 12 L is a steam sterilizer for professional instrument reprocessing in practices and clinics. It is suitable for wrapped and unwrapped instruments, hollow items and textiles, and supports Class B, S and N programs.

The SteriTrace Pro Assistant enables paperless documentation and digital batch release right on screen. Sterilization records can be saved via USB and further processed via LAN/WLAN, cloud and practice software connectivity. The included label printer also makes it easy to label sterile goods.

With a 12-litre chamber volume, 3 trays, an integrated fresh and waste water tank and an optional fixed water connection, the autoclave can be used flexibly in daily practice. This makes the Premium 12 L ideal for digitally-run dental practices, medical practices and smaller clinics focused on traceable documentation.



Scope of delivery:
Steri24 Class B Autoclave Premium 12 L
SteriTrace Pro touch panel
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
3x trays, tray holder and handle
Bacteria filter, 1x thermal paper
Versa mount for touchscreen with four mounting screws
Cold appliance power cable, mains cable, interface box, RS232 cable
2x user manuals: autoclave and SteriTrace Pro
Label printer with one roll of labels

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Premium 12 L			
Item Number	8500551-8500550	**Feed Water Tank Volume**	3 litres (for at least 5 cycles)
Chamber Volume	12 litres	**Waste Water Tank Volume**	3 litres
Chamber Dimensions	Ø 20 cm x 35 cm	**Outer Dimensions (L x W x H)**	64 x 44 x 36 cm
Sterilization Temperature	105 °C to 134 °C	**Power Connection**	230V AC, 50 Hz
Heating Power	1950 W	**Net Weight**	51 kg
Documentation	USB and SteriTrace Pro Assistant	**Connection**	RS232

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300705
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Label roll with 1000 labels	8500507
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562
Magnetic printer mount	1300160

Class B Autoclave Premium 18 L

€ 2,549

with automatic documentation

- ✓ Class B autoclave per EN 13060
- ✓ SteriTrace Pro Assistant for digital batch release
- ✓ Paperless documentation of sterilization processes
- ✓ LAN/WLAN, cloud and practice software connectivity
- ✓ Integrated fresh & waste water tank
- ✓ Label printer included
- ✓ 18-litre chamber volume with 5 trays
- ✓ Optional fixed water connection



The Steri24 Class B Autoclave Premium 18 L is a steam sterilizer for professional instrument reprocessing in practices and clinics. It is suitable for wrapped and unwrapped instruments, hollow items and textiles, and supports Class B, S and N programs.

The SteriTrace Pro Assistant enables paperless documentation and digital batch release right on screen. Sterilization records can be saved via USB and further processed via LAN/WLAN, cloud and practice software connectivity. The included label printer also makes it easy to label sterile goods.

With an 18-litre chamber volume, 5 trays, an integrated fresh and waste water tank and optional connection to a water treatment system, the autoclave is ideal for an efficient daily practice routine. This makes the Premium 18 L the right solution for dental practices, orthodontics, healthcare centres and clinics focused on networked, paperless documentation.



Scope of delivery:
Steri24 Class B Autoclave Premium 18 L
SteriTrace Pro touch panel
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
5x trays, tray holder and handle
Bacteria filter, 1x thermal paper
Versa mount for touchscreen with four mounting screws
Cold appliance power cable, mains cable, interface box, RS232 cable
2x user manuals: autoclave and SteriTrace Pro
Label printer with one roll of labels

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Pro 18 L			
Item Number	8500552-8500550	**Feed Water Tank Volume**	4 litres (for at least 5 cycles)
Chamber Volume	18 litres	**Waste Water Tank Volume**	4 litres
Chamber Dimensions	Ø 25 cm x 35 cm	**Outer Dimensions (L x W x H)**	64 x 50 x 43 cm
Sterilization Temperature	105 °C to 134 °C	**Power Connection**	230V AC, 50 Hz
Power	1950 W	**Net Weight**	53 kg
Documentation	USB and SteriTrace Pro Assistant	**Connection**	RS232

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300557
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Label roll with 1000 labels	8500507
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562
Magnetic printer mount	1300160

Class B Autoclave Premium 23 L

€ 2,799

with automatic documentation

- ✓ Class B autoclave per EN 13060
- ✓ SteriTrace Pro Assistant for digital batch release
- ✓ Paperless documentation of sterilization processes
- ✓ LAN/WLAN, cloud and practice software connectivity
- ✓ Integrated fresh & waste water tank
- ✓ Label printer included
- ✓ 23-litre chamber volume with 5 trays
- ✓ Optional fixed water connection





The Steri24 Class B Autoclave Premium 23 L is a steam sterilizer for professional instrument reprocessing in practices, healthcare centres and clinics. It is suitable for wrapped and unwrapped instruments, hollow items and textiles, and supports Class B, S and N programs.

The SteriTrace Pro Assistant enables paperless documentation and digital batch release right on screen. Sterilization records can be saved via USB and further processed via LAN/WLAN, cloud and practice software connectivity. The included label printer also makes it easy to label sterile goods.

With a 23-litre chamber volume, 5 trays, an integrated fresh and waste water tank and an optional fixed water connection, the autoclave is ideal for a high volume of instruments. This makes the Premium 23 L the powerful solution for clinics, surgical centres, large dental practices and healthcare centres focused on digital batch release.

Scope of delivery:
Steri24 Class B Autoclave Premium 23 L
SteriTrace Pro touch panel
2x quick-connect adapters and 1x hose for manually draining the feed and waste water tanks
5x trays, tray holder and handle
Bacteria filter, 1x thermal paper
Versa mount for touchscreen with four mounting screws
Cold appliance power cable, mains cable, interface box, RS232 cable
2x user manuals: autoclave and SteriTrace Pro
Label printer and one roll of labels

Click here for our current offers:

Specifications

Steri24 Class B Autoclave Premium 23 L			
Item Number	8500553-8500550	**Feed Water Tank Volume**	4 litres (for at least 5 cycles)
Chamber Volume	23 litres	**Waste Water Tank Volume**	4 litres
Chamber Dimensions	Ø 25 cm x 45 cm	**Outer Dimensions (L x W x H)**	68 x 50 x 43 cm
Sterilization Temperature	105 °C to 134 °C	**Power Connection**	230V AC, 50 Hz
Power	1950 W	**Net Weight**	55 kg
Documentation	USB and SteriTrace Pro Assistant	**Connection**	RS23

Recommended Accessories

Item	Item Number
Bacteria filter (1 pc.)	1300705
Helix test (with 100 test strips)	8500434
Steam Plus sterilization test (with 100 test strips)	8500436
Label roll with 1000 labels	8500507
Fixed water connection	8500564
AQUAdem 40 – Water Treatment System	8500562
Label printer holder	1300160

Wide Range of Accessories

For an optimal workflow

Tests & Steam Indicators

For use in all Steri24 Class B Pro autoclaves

Accessories	Item Number	Price
Helix test (incl. 100 test strips) - suitable for hollow items	8500434	€ 45.00
Bowie & Dick test (20 test strips) - suitable for textiles	8500438	€ 65.00
Steam Plus sterilization test (100 test strips) - not suitable for hollow items	8500436	€ 54.00


8500434


8500438


8500436

Paper

The perfect choice for use in the built-in printer of Steri24 Class B Pro autoclaves. Each roll provides room for about 50 protocols.

Accessories	Item Number	Price
Thermal paper rolls (10 pcs.)	1300742	€ 12.50


1300472

Filter

For the Steri24 Class B Pro autoclave. The filter must be replaced every 150 cycles. With fewer than 150 cycles per year, the filter must be replaced once a year.

Accessories	Item Number	Price
Bacteria filter (1 pc.)	1300557	€ 19.90



1300557

Inserts and Base Plates

As an option to the included tray holder and integrated trays, we also offer aluminium trays, an insert for wrapped instruments, and base plates for use with sterilization boxes.

Description	Item Number	Price
Aluminium Tray 12 L	1300907	€ 39.00
Aluminium Tray 18 L	1300908	€ 44.00
Aluminium Tray 23 L	1300909	€ 49.00
Insert for wrapped instruments 23 L	8500548	€ 35.00


8500548


1300907



Water Supply

High-quality, clean, demineralized water is especially important for autoclaving. Poor water quality leads to deposits in the autoclave chamber. That is why we offer two water treatment systems for all our autoclaves.

Water distiller

This lets you save money compared to buying demineralized water.

Accessories	Item Number	Price
Water distiller, 1.5 litres / hour	8500482	€ 99.00
Activated carbon filter for water distiller (12 pcs.)	8500482AK	€ 51.00
AQUAdem 40 – Water Treatment System	8500562	€ 299.00
Spray gun for AQUAdem 40	1300565	€ 49.00
Mixed-bed resin cartridge for AQUAdem 40 (1 pc.)	1300553	€ 39.00
Fixed water connection at the autoclave	8500564	€ 499.00



8500482


1300565



8500562

Accessories for Class B Autoclave Premium series

To complement our Class B Premium series, we offer a label printer for labelling your sterile goods. The barcode scanner allows batch numbers to be scanned in.

Label Printer and Barcode Scanner

For comprehensive traceability of your sterile goods.

Accessories	Item Number	Price
Label printer	8500502	€ 224.00
Label roll with 1000 labels	8500507	€ 25.00
Label printer holder incl. magnets	1300160	€ 49.00
Barcode scanner	8500503	€ 121.00


8500502

8500507


1300160


8500503

CLEANING
ULTRASONIC CLEANERS



AI-generated image

Steri24 Ultrasonic Cleaners

from € 330

Professional ultrasonic cleaning for the highest standards

Steri24 ultrasonic cleaners stand for reliable, efficient and thorough cleaning in professional work environments. Designed for daily use in practices, clinics, labs, workshops and technical operations, they support effective pre-cleaning of instruments, tools, components and delicate small parts.

With an ultrasonic frequency of 40,000 Hz, the units loosen even stubborn dirt, deposits and residue in hard-to-reach places. This enables especially even, material-friendly cleaning – ideal for applications where precision, hygiene and reliable results are essential.

The Steri24 range offers various tank sizes from 3.2 to 22 litres, covering a broad range of applications: from compact units for smaller cleaning tasks to large-volume solutions for demanding cleaning needs and larger instruments or workpieces.

The adjustable heating from 20 °C to 80 °C and the mechanical timer from 1 to 30 minutes allow flexible adaptation to different cleaning processes. The included stainless steel basket ensures safe placement of the items to be cleaned and simplifies handling in professional daily use.

Thanks to their robust construction, simple operation and versatile applications, Steri24 ultrasonic cleaners are a reliable solution for all areas where thorough cleaning, efficiency and a smooth workflow are required.





Your benefits at a glance

- ✓ Professional ultrasonic cleaning at 40,000 Hz
- ✓ Various tank sizes for individual requirements
- ✓ Adjustable heating from 20 °C to 80 °C
- ✓ Mechanical timer from 1 to 30 minutes
- ✓ Stainless steel basket included
- ✓ Robust construction for daily use
- ✓ Suitable for instruments, tools, components, jewelry, glasses and more
- ✓ Ideal for practices, clinics, labs, workshops and technical applications

Steri24 ultrasonic cleaners – for clean results, efficient workflows and professional cleaning processes.

Click here for our current offers:

Specifications


Tank Capacity	Tank Size (L x W x H)	Ultrasonic Power	Item Number	Price
3.2 L	240 x 135 x 100 mm	120 W	8500537	€ 330.00
6.5 L	300 x 150 x 150 mm	180 W	8500538	€ 570.00
10 L	300 x 240 x 150 mm	240 W	8500539	€ 747.00
15 L	330 x 300 x 150 mm	360 W	8500540	€ 990.00
22 L	500 x 300 x 150 mm	480 W	8500541	€ 1,227.00

Recommended Accessories

Item	Item Number
Multisteril CD Ctrl Alt Solution for ultrasonic unit, 1 litre	8500458 € 45.00



DISINFECTION
THERMAL DISINFECTOR



AI-generated image

Miele Thermal Disinfector

PWD 8534 countertop – disinfection & drying

from € 7,887



- ✓ Thermal disinfection per EN ISO 15883
- ✓ DryPlus hot-air drying with HEPA-14 filter
- ✓ Two integrated dosing pumps for automatic dosing
- ✓ Full-glass door for visual inspection of the wash chamber
- ✓ 8 programmable user programs with direct selection buttons
- ✓ Integrated water softener and triple water filter system
- ✓ Capacity: 8 hollow instruments plus additional instruments

The Miele PWD 8534 countertop thermal disinfector is optimized for use in small practices and offers a state-of-the-art solution for efficient reprocessing of instruments with or without hollow items.

DryPlus hot-air drying with HEPA-14 filter ensures thorough drying, while two integrated dosing pumps allow automatic dosing of cleaning and care agents. The full-glass door allows visual inspection of the wash chamber at any time.

With a 73-litre wash chamber volume, 8 programmable user programs with direct selection buttons, an integrated water softener and a triple water filter system, the PWD 8534 offers a compact solution for small practices. Its dimensions are 555 x 645 x 575 mm (H x W x D) with a weight of approx. 55 kg.



Scope of delivery:
Miele PWD 8534 countertop thermal disinfector with full-glass door, Water softener, basket accessories selectable separately

Click here for our current offers:

Specifications

Miele Thermal Disinfector PWD 8534			
Item Number	PWD8534	**Dimensions (H x W x D)**	555 x 645 x 575 mm
Wash Chamber Volume	73 litres	**Weight**	approx. 55 kg
Heating Power	2.6 kW	**Programs**	8, directly selectable

Recommended Accessories

Item	Price
APWD 060 base basket	€ 150
APWD 079 injector basket	€ 890
E 491 insert for rotary instruments	€ 595
E 521/2 insert for 7 forceps	€ 154
Validation (DE/AT/CH)	on request
Commissioning with instruction	on request

Miele Thermal Disinfector

PWD 8682 – disinfection & drying

from € 8,628



- ✓ Thermal disinfection per EN ISO 15883
- ✓ Final rinse temperature up to 93 °C
- ✓ EcoDry drying with AutoOpen function
- ✓ Laser-welded wash chamber, 145-litre volume
- ✓ M Touch Flex touchscreen with 8 freely programmable programs
- ✓ Ethernet/WiFi networking ex works
- ✓ Capacity: 44 handpieces or suction cannulas per load

The Miele PWD 8682 is a powerful thermal disinfector for professional instrument reprocessing in dental practices. Thermal disinfection is carried out with hot water and demineralized water connections in a validatable process per EN ISO 15883.

Integrated wash pressure and spray arm monitoring immediately detects pressure drops or blockages, while the variable-speed heater pump optimizes wash pressure for each program block. EcoDry with automatic door opening (AutoOpen) ensures faster, thorough drying.

With a 145-litre wash chamber volume, 8 freely programmable programs, a color M Touch Flex touchscreen and factory-fitted Ethernet/WiFi networking, the PWD 8682 holds e.g. 44 handpieces or 44 suction cannulas per load. Its dimensions are 835 x 598 x 598 mm (H x W x D).



Scope of delivery:
Miele PWD 8682 thermal disinfector with touchscreen and WiFi networking, Baskets selectable separately

Click here for our current offers:

Specifications

Miele Thermal Disinfector PWD 8682			
Item Number	8500582-LW	**Dimensions (H x W x D)**	835 x 598 x 598 mm
Wash Chamber Volume	145 litres	**Shortest program time**	40 min
Max. final rinse temperature	93 °C	**Programs**	8, freely programmable

Recommended Accessories

Item	Price
Dental basket set APWD 212	€ 2,837
Miele liquid detergents APWd 212 Dental	€ 62
Validation package incl. commissioning (DE/AT)	€ 1,549
Commissioning incl. travel (DE/AT)	€ 699
Upper/lower basket range (A 101–A 105/1)	on request
Trolley range (A 201–A 208)	on request

Matching accessories

For your practice

Detergents

the right detergents for your unit

Accessories	Capacity	Item Number	Price
Pro Care Dent 10 MA	5 litres	8500584	€ 62.00
Pro Care Dent 30 C	5 litres	8500585	€ 71.00
Pro Care Dent 40	1 litre	8500586	€ .00
Pro Care regeneration salt 61	6 kg	8500572	€ 13.00










Miele basket set APWD 212 Dental


Included accessories	8500583	€ 2,837.00
A105/1 upper basket		
E473/2 mesh tray with lid (for various small parts)		
E337/1 insert for upright instruments		
A836 injector nozzle set		
A803 holder for 3 handpieces		
A800 filter tube		
A804 cleaning brush		
A151 lower basket		


A800


A105/1


A804


E337/1


E337/1


A151


A836


A803

Every application is different. That is why we put together your accessory package individually and will gladly prepare a suitable quote for you.

PACKAGING
SEALING DEVICES, POUCHES & BOXES



AI-generated image

CertoSeal Pro Touch

€ 475

Sealer with touchscreen

- ✓ Making custom pouches from roll stock
- ✓ Integrated cutting function
- ✓ Adjustable temperature from 90 °C to 210 °C
- ✓ Roll holder for rolls up to 30 cm wide
- ✓ Counter for sealed packages



The CertoSeal Pro Touch is a flexible, affordable entry-level sealer for professional instrument reprocessing. Existing pouches can be securely sealed, or custom packaging can be made to size from roll stock.

The integrated cutting function and adjustable temperature make the unit suitable for different packaging materials and everyday practice requirements. Simple operation and a counter for sealed packages support an efficient, well-organized workflow.

The roll holder is suitable for rolls up to 30 cm wide. This makes the CertoSeal Pro Touch ideal for general practices, veterinary practices, medical, eye and hair clinics, and facilities without a validation requirement.

Scope of delivery:
CertoSeal Pro Touch sealer
Replacement blade
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

	CertoSeal Pro Touch
Product Model	CertoSeal Pro Touch
Dimensions (L x W x H)	28.3 x 51.7 x 20.5 cm
Temperature	90 °C to 210 °C
Power Connection	230V AC, 50 Hz

Recommended Accessories

Item	Item Number
Sterilization paper roll 7.5 cm x 200 m	8500456
Sterilization paper roll 10 cm x 200 m	8500455
Sterilization paper roll 15 cm x 200 m	8500454
Sterilization paper roll 20 cm x 200 m	8500453
Sterilization paper roll 25 cm x 200 m	8500452

CertoSeal 200

€ 999

Validatable continuous sealer with cutter

- ✓Validatable & documentable
- ✓Automatic pouch feed
- ✓Manual cutting function for roll stock
- ✓7-inch touchscreen
- ✓USB port for sealing process documentation



The CertoSeal 200 is a validatable continuous sealer for professional instrument reprocessing in practices and clinics. Pouches are fed in automatically and sealed securely with consistent, reproducible quality.

The large 7-inch touchscreen makes the unit simple and intuitive to operate. The integrated USB interface enables reliable documentation of sealing processes, supporting the requirements for traceable reprocessing.

Thanks to the roll holder for rolls up to 30 cm wide and the manual cutting function, packaging can be prepared flexibly and cleanly. This makes the Certo-Seal 200 the cost-conscious, validatable solution for dental practices, medical practices, veterinary practices, eye and hair clinics, and other medical facilities.



Scope of delivery:
CertoSeal 200 sealer
Replacement blade
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

	CertoSeal 200
Product Model	CertoSeal 200
Dimensions (L x W x H)	54 × 23.5 × 21.5 cm
Power Connection	230V AC, 50 Hz
Sealing Temperature	Adjustable from room temperature to 220°C
Seal Seam Width	12mm
Heating Speed	Rises to 180°C in 40 seconds
Margin	0–40 mm
Power	400W
Net Weight	16 kg

Recommended Accessories

Item	Item Number
Sterilization paper roll 7.5 cm x 200 m	8500456
Sterilization paper roll 10 cm x 200 m	8500455
Sterilization paper roll 15 cm x 200 m	8500454
Sterilization paper roll 20 cm x 200 m	8500453
Sterilization paper roll 25 cm x 200 m	8500452

CertoSeal 300

Validatable continuous sealer with **autocutter**

€ 1,299

- ✓ Validatable & documentable
- ✓ Automatic pouch feed
- ✓ Manual cutting function for roll stock
- ✓ 7-inch touchscreen
- ✓ USB port for sealing process documentation





The CertoSeal 300 was developed for facilities that place particular value on a fast and convenient packaging process during instrument reprocessing. The validatable continuous sealer feeds in pouches automatically, seals them reliably and automatically cuts the packaging material to the required length.

Operation is convenient via a clear 7-inch touchscreen. All relevant functions can be selected directly, while the integrated USB interface enables complete documentation of sealing operations. This allows the CertoSeal 300 to support a transparent, standards-compliant reprocessing process.

The combination of automatic feed, automatic cutting function and roll holder for packaging rolls up to 30 cm wide significantly reduces repetitive work steps. This saves time, simplifies handling and ensures a smooth workflow in daily practice and clinic operations.

This makes the CertoSeal 300 the convenience-focused solution for dental practices, medical practices, veterinary practices, eye and hair clinics, and other medical facilities with regular packaging needs.

Scope of delivery:
CertoSeal 300 sealer
Replacement blade
Cold appliance power cable
User manual

Click here for our current offers:

Specifications

	CertoSeal 300
Product Model	CertoSeal 300
Dimensions (L x W x H)	54 × 23.5 × 21.5 cm
Power Connection	230V AC, 50 Hz
Sealing Temperature	Adjustable from room temperature to 220°C
Seal Seam Width	12mm
Heating Speed	Rises to 180°C in 40 seconds
Margin	0–40 mm
Power	400W
Net Weight	16 kg

Recommended Accessories

Item	Item Number
Sterilization paper roll 7.5 cm x 200 m	8500456
Sterilization paper roll 10 cm x 200 m	8500455
Sterilization paper roll 15 cm x 200 m	8500454
Sterilization paper roll 20 cm x 200 m	8500453
Sterilization paper roll 25 cm x 200 m	8500452

Wide Range of Accessories

For an optimal workflow



Sterilization Rolls

Sterilization paper for use with sealing devices.


Accessories	Item Number	Price
Sterilization paper roll, 7.5 cm wide, 200 m long	8500456	€ 18.00
Sterilization paper roll, 10 cm wide, 200 m long	8500455	€ 22.00
Sterilization paper roll, 15 cm wide, 200 m long	8500454	€ 33.00
Sterilization paper roll, 20 cm wide, 200 m long	8500453	€ 45.00
Sterilization paper roll, 25 cm wide, 200 m long	8500452	€ 56.00

Sterilization Paper Pouches

Self-sealing sterilization paper pouches, also known as sterilization pouches or laminated pouches. For use without a sealing device.

Accessories	Item Number	Price
Self-sealing pouches 90 x 250 mm (200 pcs.)	8500449	€ 13.00
Self-sealing pouches 140 x 260 mm (200 pcs.)	8500448	€ 12.00
Self-sealing pouches 190 x 330 mm (200 pcs.)	8500447	€ 20.00
Self-sealing pouches 300 x 450 mm (200 pcs.)	8500446	€ 65.00





The benefits of Steri24:

Free delivery in DE, AT and CH
for orders of € 1,200 net or more

2–3 year warranty
2 years' warranty as standard; 3 years with a maintenance contract

24-hour emergency hotline
+49 545 263 499 31 - Steri24 service location near you

On-site service
Commissioning, maintenance and validation of your equipment on-site by our service partners

Base plate for sterilization containers

Accessories	Item Number	Price
Base Plate 12 L	1300910	€ 39.00
Base Plate 18 L	1300911	€ 44.00
Base Plate 23 L	1300912	€ 49.00

Sterilization Containers

Aluminium sterilization containers for safe sterilization and storage.

Accessories	Size (L x W x H)	Item Number	Price
Aluminium Sterilization Container Mini S	300 x 140 x 40 mm	1300923	€ 99.00
Aluminium Sterilization Container Mini M	300 x 140 x 70 mm	1300924	€ 119.00
Aluminium Sterilization Container Mini L	300 x 140 x 100 mm	1300925	€ 129.00
Aluminium Sterilization Container Dental S	310 x 190 x 40 mm	1300926	€ 119.00
Aluminium Sterilization Container Dental M	310 x 190 x 65 mm	1300927	€ 129.00
Aluminium Sterilization Container Dental L	310 x 190 x 130 mm	1300928	€ 149.00

Accessories for Sterilization Containers


Single-use or reusable filters to fit our Aluminium Mini and Dental sterilization containers.

Accessories	Item Number	Price
Single-use filters (250 pcs.)	1300930	€ 15.00
Teflon reusable filters (2 pcs.)	1300931	€ 39.00
Blue indicator seal (1,000 pcs.)	1300952	€ 159.00









Your strong partner!

Steri24 is known not only for its autoclaves, but also for its service and additional offerings. We advise you on sterilization processes, help you find the right autoclave, and remain at your side after purchase.
On request, we can carry out complete qualification and validation for individual Steri24 models.

Services

On-site equipment maintenance	For all Steri24 units	Item number: 1250006

Under legal requirements, your autoclave, thermal disinfector and sealer should be serviced every 24 months (for autoclaves, also after 1000 sterilization cycles at the latest).

You can easily order maintenance at a fixed price in our online shop.
Tip: sign a maintenance contract with us. We will make sure your autoclave is always serviced on time.
This offer, including travel, is available in Germany, Austria, Switzerland, Luxembourg and Liechtenstein. Our service package covers the cost of maintenance, travel and standard replacement parts (door seal, bacteria filter, safety valve). You will only be invoiced once maintenance is complete.

On-site validation	For all Steri24 autoclaves & disinfectors	Item number 1250003

Validation (installation, operational and performance qualification) and revalidation (repeat performance qualification) are carried out based on the relevant parts of DIN EN ISO 15883 and the current guidelines from DGKH, DGSV and RKI for the validation and routine monitoring of mechanical cleaning and disinfection processes, as well as the guidelines for steam sterilization and seal seam processes.
The current recommendations of the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM), “Requirements for Hygiene in the Reprocessing of Medical Devices”, are also taken into account.

Tip: sign a validation contract with us. We will make sure your autoclave is always revalidated and stays ready for use.
With a validation contract, you also get a price advantage thanks to better planning for on-site validation.

Repair service	For all Steri24 units	Item number: REPAT

You get a standard two-year warranty with us. During this time, we carry out repairs free of charge under the warranty terms.

If the worst happens, you can easily order a repair at the base price in our online shop.
This offer, including travel, is available in Germany, Austria, Switzerland, Luxembourg and Liechtenstein.
Our repair package covers the cost of the repair, travel and standard replacement parts (door seal, bacteria filter, safety valve).



AI-generated image

The benefits of Steri24:

Free delivery in DE, AT and CH
for orders of € 1,200 net or more

24-hour emergency hotline
+49 545 263 499 31 – Steri24 service location near you

2–3 year warranty
2 years’ warranty as standard; 3 years with a maintenance contract

On-site service
Commissioning, maintenance and validation of your equipment on-site by our service partners

Questions? We are always happy to help!

For more information, visit www.steri24.com or contact our office directly.

Contact Details

Address
Steri24 – a brand of
CertoClav Sterilizer GmbH
Peintner Straße 10
A-4060 Leonding

Phone
+43 732 674 278 22

Email
support@certoclav.com

Company Register
Regional Court Linz
FN122912d

VAT No.
ATU22821702

CEO
Michael Simon

For advice, contact us at
support@certoclav.com
Tel. +43 732 674 278 22
Monday – Thursday: 8:00 am – 4:00 pm
Friday: 8:00 am – 1:00 pm

Order around the clock at our Steri24 online shop:
www.steri24.com

Notes

The warranty terms are described in the respective certificates.
Deliveries are made under our terms and conditions, available at: https://www.steri24.com.
Product images may differ slightly from the product. Images are not to scale relative to each other.
All prices are net plus statutory VAT.
Actual prices may vary due to customary industry promotions and variable purchase prices.

Imprint

© Certoclav Sterilizer GmbH
The images, information and technical data in this publication reflect the status at the time of printing. The author assumes no liability for the timeliness, accuracy, completeness or quality of the information provided.

Notes:

Orders:

by email: support@certoclav.com
by phone: +43 732 674 278 22
in our online shop: www.steri24.com



AI-generated image

Steri24

Order now
+43 732 674 278 22


By email
support@certoclav.com

online shop
www.steri24.de

ClavioSoft

- Fully automatic documentation
- Over 200 devices connectable
- RKI-compliant complete Documentation
- Optimized for hygiene audits

For more information, visit www.steri24.de



Address
Steri24 - a brand of
CertoClav Sterilizer GmbH
Peintner Straße 10
A-4060 Leonding

Phone
+43 732 674 278 22

Email
support@certoclav.com

All technical data apply exclusively to unloaded units in standard configuration at an ambient temperature of 22 ± 3 °C and a mains voltage fluctuation of ±10 %. Subject to technical changes.

www.steri24.com

01.2027